12 Divers

Various

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Benefit of an unguual prosthesis and setting technique.
Dr S. Guero

1. In the Emergency Room

Avulsion of the unguual tablet is accompanied by an impairing of the nail bed with two unsatisfactory developments: hyper-budding and keratinization.


Hyper-budding

The nail bed being very vascularized, botryomycomas are frequent. They are sometimes resistant to all medical treatments (corticosteroids, nitrating) and, in some cases, surgical excision is necessary. Replacing the tablet or placing an unguual prosthesis mechanically prevents hyper-budding.

Keratinization

Henceforth, we know that the nail bed is not a sterile matrix (Fig. 1). According to the individual, keratin production varies between one and ten percent of the thickness of the unguual tablet. In the case of an unguual avulsion, after the scarring-over phase, the nail bed reacts to air exposure by starting to produce keratin.

In the best cases, the newly-formed nail is supple and butts upon the hardened nail bed. The nail slows its progress from proximal to distal and forms a cushion proximal to the nail bed. This cushion will only then thicken very slowly. The unguual tablet will arrive at the hyponychium in six months, instead of three or four, and will be very unattractive. Six more months must pass before it becomes refined.

2. "Cold" surgery for unguual dystrophies

Consultation requests for unguual dystrophies are very frequent, since the esthetics of the hands and feet are important to patients. Of traumatic, infectious, or congenital origin, unguual dystrophies may be improved, even healed, provided that the regrowth of the new nail is guided in order to avoid physiological mishaps such as the ones mentioned above. Because this removal is painful, prescribe numerous bandages and sometimes cause unguual dystrophies at the given site. The remedy of an unguual prosthesis, whose make-up is close to the nail morphology, simplifies considerably the method of this surgery. The unguual prostheses are thus indispensable in repairing surgery :

- in urgent cases where the extracted nail is no longer available,
- "cold", in order to avoid the removal of a healthy nail.

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Technical note

The ungual prosthesis is slightly larger than normal in order to enable resizing and adaptation to all nail forms. Its concavity reproduces the form of a normal nail. The fundamental quality is that its underlying side is striated longitudinally like a natural nail. The nail must be attached laterally by two separate points, or with a U-point (Fig. 2a and Fig. 2b).

On the other hand, it is not advised to place a point in the dorsal skin at the risk of arriving in the matrix and leaving a permanent longitudinal ridge (Fig. 2c). This setting technique advocated by certain authors must be done away with.


Moreover, the prosthesis will only be inserted under the skin during "cold" surgery, in reconstructing the groove of the matrix.

In urgent cases, we suggest inserting a gap between the dorsal skin and the sub-cutaneous part of the nail (soft nail), which lessens the risks of peringual (?) whitlow. Even if a prosthesis delivered sterile is better tolerated than an avulsed nail, it is imperative to insure proper drainage of the hæmatoma.

Bibliographie / Bibliography

- **Achten G.** : l'ongle normal, J. Med. Esth. Chir. Derm., XV: 1988; XV: 193-200.
- **Baran R.** : in l'ongle, Monographie du GEM N°9. Expansion Scientifique Française, 1978; Paris
- **Baran R., Bureau H.** : Congenital malalignment of the big toe-nail as a cause of ingrowing toe-nail in infancy. Pathology and treatment (a study of thirty cases). Clinical and Experimental Dermatology, 1983; 4: 309-313.
- **Baran R., Dawber R.P.R.** : Diseases of the nails, blackwell Scientific Publications, 1984; Oxford.
- **Guero S., Gulchard S., Fraitag S.** : Ligamentary structure of the base of the nail. Surgical and Radiological Anatomy, 1984; 16: 47-52.
- **Johnson M., Comaish J. S., Shuster S.** : Nail is produced by the normal nail bed : a controversy resolved. Br. J. Dermatol., 1991; 125: 27-29.
- **Lewis B.L.** : (1954) Microscopic studies of fetal and mature nail surrounding soft tissue. Arch. Dermatol., 1954; 70: 732-747.
- **Shepard G. H.** (1954): Treatment of nail bed avulsions with split thickness nail bed grafts. J. Hand Surg., 1983; 8: 45-54.
- **Zalas N.** : The Nail, Jamaica. NY, Spectrum Publishers, 1980.
- **Zook E. G., Van Beek A. L., Russell R. C. et al.** : Anatomy and physiology of the perionychium. A review of the litterature and anatomic study. J. Hand Surg., 1986; 5: 528-536.
- **Zook E. G.** : The perionychium. Anatomy, physiology and care of injuries. Clin. Plast. Surg., 1984; 8: 21-31.
- **Zook E. G., Guy R. J., Russell R. C.** : A study of nail bed injuries: Causes, treatment and prognosis. J. Hand Surg., 1984; 9A: 247-252.
- **Zook E. G.** : The perionychium. Hands Clinics. 1990; 6: N° 1.

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Cupula prostheses are interposition prostheses between the metatarsal head and the re-cut phalanx, set to the capsule with catgut or 000 nylon stitches.

These prostheses do not have an intra-phalangeal fixing eye-bolt, thus causing no bone reactions, such as osteolysis. Being thin (8/10 millimeters), they cause no intolerance or inflammatory reactions.

Acting as a meniscus, they prevent bone ankylosis and insure the continuation of articular amplitude.

They are especially appropriate in surgical therapy of arthroses like Hallus Rigidus, Freiberg Koehler, and in iatrogenic and post-traumatic arthroses.

Interphalangeal use of the Regnault Cupula

In the case of a forty-year-old man: painful stiffening, unaligned proximal inter-phalanx of the second toe after traumatic dislocation, reduced in emergency, unimmobilized.

- **Operation at 4 months**: arthroplasty of the proximal inter-phalanx with stringy interposition.
 - **On the fourth day**, direct shock to the second toe, the results are very painful with toe distension, progressive stiffening, appearance of capsular ossification, pain caused by the least mobilization of the second toe.
 - **Six months after the operation**: painful stiffness of the proximal inter-phalanx. The patient wants to regain mobility; to avoid the impaction of P2 on P1, the interposition of a Regnault cupula re-cut to the dimensions proximal inter-phalanx is planned, after the recut of P1 and P2, and to exeresse the capsular ossifications.
- Operation realised in september 1990.*
- **Six months after the second operation**, the second toe has recovered its normal volume, exhibits no pain, flexion is satisfactory, extension is normal. Only a few vascular troubles persist. The cupula will be removed only if it poses a problem.

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Ø	Référence
10 mm	61.001.10
14 mm	61.001.14
18 mm	61.001.18
22 mm	61.001.22

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Indication

The titanium button is used for trans-osseus sutures of the caul to support the node when the bone is soft. It may be used for all trans-osseous sutures by U-threading in which the bone is soft (fractures of the trocheter, tuberosity, etc.).

Thread to use

Any nylon thread, but it must be sufficiently resistant so as not to be quickly cut in two by the button.

Why titanium ?

In order to permit modern post-operative imaging via a scanner or M.R.I.

Noteworthy: The width of the button for the shoulder is very small, so there is no conflict with the trochiter and no protrusion.

Forage Resection Osseuse Percutanees (FROP) Percutaneous Bone Resection Drilling (PEBORD)

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Contrôles: J 45, 6 mois, 1 an

For a long time, the treatment of osteoid osteoma involved a "block" osseous excision that was damaging and disproportionate with respect to the size of the lesion.

The difficulties in locating the nidus during surgery and the desire to limit the bone resection as much as possible contributed to the development of the percutaneous excision technique with scanographic location and guidance.

Advantages of FROP

(percutaneous bone resection drilling)

Precision of the location

- Very limited nature of the excision (not very invasive procedure)
- Little risk of post-operative fracture
- Short-term hospitalization (48 to 72 hours)
- Rapid recovery of ability to use legs (with canes)

Resumption of normal activity is possible about 2 months after the treatment

Percutaneous resection technique: FROP (percutaneous bone resection drilling)

General anesthesia

Collaboration of surgeon and radiologist

Duration of the procedure: 45 to 90 minutes (locating and resection)

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1 Accurate locating of the nidus




The examination is done under local or general anesthesia. The lesion is located by making meeting tomographic cuts in order to define the ballistic strategy. A cut centered on the nidus is made to estimate its diameter and depth. In order to facilitate the drilling, the member is positioned so that the angle of attack is as direct as possible (*as vertical as possible*).

2 Percutaneous resection


A two centimeter-long cutaneous incision is made at the chosen point of entry marked on the skin. The following steps are then carried out:

- 1) The trocar with its mandrel is placed facing the lesion (*image 1*)
- 2) The mandrel is replaced by a pin, which must go into the nidus. The correct position of the pin is verified by several tomographic cuts. (*image 2*)
- 3) Once the pin has been properly positioned facing the nidus, the trocar is introduced up to the surface of the bone. A cover is slipped onto the trocar, which is then removed. (*image 3 and 4*)
- 4) A 6 mm diameter drill bit is introduced on the pin. The drilling is done with a motorized unit. The same action is then done with a 9 mm diameter drill bit. (*image 5*)
- 5) After the drilling, the zone of the nidus is approached by the trephine, which extracts the bone sample containing it. (*image 6 and 7*)
Then milling is necessary to destroy the region adjacent to the nidus to prevent a possible relapse.
- 6) After the procedure, a scanographic check is done to verify complete excision and the absence of complications. If the excision is incomplete, a new drilling is done.

Few severe complications have been reported and the technique is less expensive than surgical resection.

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Taraudage du trou percé au préalable
Tapping of the hole bored as a preliminary

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The trocart allows the drilling & the injection of osseous marrow.

The trocart allows a percutaneous use with an entrance point in the area trochantérienne, an aiming of the femoral head under amplifier of brightness, and an intraosseous penetration with the hammer.

When the trocart is placed, one can leave in place the sleeve and withdraw the chuck. This sleeve then allows the injection of osseous marrow or another liquid substance using a syringe of 10 mL.

Then: reinstallation of the chuck and ablation.

This device allows the taking away and the injection of marrow in same operational time.

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Instrument for making punctures composed of:

- a metallic socket
- a mandrel: metallic rod that has a triangular tip with sharp edges at its end that goes 10 mm beyond the socket.

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